



ACCESSORIES

BOLUS – GEL, SOLID +PELLETS

GEL BOLUS

Gel Bolus products are flexible, tissue equivalent build-up materials which, when applied on a patient, increase the surface dose uniformly.

Gel Bolus is available with a durable and cleanable layer of thin plastic skin or as a skinless type, which is sticky for conformity to body contours and easily cut to specialized shapes as needed.



SOLID BOLUS

Solid Bolus is ideal for irregularly-surfaced anatomical areas that will benefit from rigid conformal bolus material.

BOLUS PELLETS

Bolus pellets, once softened in warm water or dry heat, can be shaped and bonded to almost any required area to form a permanent fixture.

Not available in all markets. Please ask your sales representative.

Thickness	Size	With Skin	Skinless
0.3cm	15cm x 15cm	g-BOLUS-A-15x15x0.3	
	30cm x 30cm	g-BOLUS-A-30x30x0.3	
	30cm x 56cm	g-BOLUS-A-30x56x0.3	
	40cm x 40cm	g-BOLUS-A-40x40x0.3	
0.5cm	15cm x 15cm	g-BOLUS-A-15x15x0.5	
	20cm x 20cm	g-BOLUS-A-20x20x0.5	
	30cm x 30cm	g-BOLUS-A-30x30x0.5	g-BOLUS-B-30x30x0.5
	30cm x 40cm	g-BOLUS-A-30x40x0.5	
	30cm x 45cm	g-BOLUS-A-30x45x0.5	
	30cm x 50cm	g-BOLUS-A-30x50x0.5	
	30cm x 56cm	g-BOLUS-A-30x56x0.5	
	40cm x 30cm		g-BOLUS-B-40x30x0.5
	40cm x 40cm	g-BOLUS-A-40x40x0.5	g-BOLUS-B-40x40x0.5
	40cm x 56cm	g-BOLUS-A-40x56x0.5	
1.0cm	50cm x 50cm	g-BOLUS-A-50x50x0.5	g-BOLUS-B-50x50x0.5
	10cm x 15cm	g-BOLUS-A-10x15x1.0	
	15cm x 15cm	g-BOLUS-A-15x15x1.0	
	20cm x 20cm	g-BOLUS-A-20x20x1.0	
	30cm x 30cm	g-BOLUS-A-30x30x1.0	g-BOLUS-B-30x30x1.0
	30cm x 50cm	g-BOLUS-A-30x50x1.0	
	30cm x 56cm	g-BOLUS-A-30x56x1.0	
	40cm x 40cm	g-BOLUS-A-40x40x1.0	g-BOLUS-B-40x40x1.0
	40cm x 56cm	g-BOLUS-A-40x56x1.0	g-BOLUS-B-40x56x1.0
	50cm x 50cm	g-BOLUS-A-50x50x1.0	
1.5cm	56cm x 56cm	g-BOLUS-A-56x56x1.0	g-BOLUS-B-50x50x1.0
	15cm x 15cm	g-BOLUS-A-15x15x1.5	
2.0cm	30cm x 30cm	g-BOLUS-A-30x30x1.5	g-BOLUS-B-30x30x1.5
	10cm x 15cm	g-BOLUS-A-10x15x2.0	
3.0cm	15cm x 15cm	g-BOLUS-A-15x15x2.0	
	30cm x 30cm	g-BOLUS-A-30x30x2.0	g-BOLUS-B-30x30x2.0
	30cm x 30cm	g-BOLUS-A-30x30x3.0	g-BOLUS-B-30x30x3.0
4.0cm	30cm x 30cm		g-BOLUS-B-30x30x4.0

Thickness	Size	Solid Bolus
4.8mm	20cm x 23cm	S-BOLUS-4.8x20x23
	30cm x 45cm	S-BOLUS-4.8x30x45
	43cm x 43cm	S-BOLUS-4.8x43x43

453.6 gram (1lb) bottle of pellets	S-BOLUS-P1
------------------------------------	-------------------