



FREEDOMX™

ACCESSORIES

PNEUMATIC ABDOMINAL INDEXABLE COMPRESSION BELT SYSTEM

The FreedomX™ Indexable Compression Belt is an easy to use pneumatic abdominal compression system for managing movement and tumor motion while maintaining maximum comfort to the patient. The belt indexes to the FreedomX™ Overlay System for clear, easy accuracy.



COMP-BLT-SX	FreedomX™ Indexable Small Compression Belt System with 1 Sphygmomanometer
COMP-BLT-LX	FreedomX™ Indexable Large Compression Belt System with 1 Sphygmomanometer
COMP-KIT-AX	FreedomX™ Indexing Compression Belt Kit Includes: Small and Large Belts with 1 Sphygmomanometer

STANDARD ABDOMINAL COMPRESSION BELTS

The FreedomX™ Compression Belt is an easy to use pneumatic abdominal compression system for managing movement and tumor motion while maintaining maximum comfort to the patient. The belt can be used in conjunction with the SBRT modules of the FreedomX™ System, SaBella Flex™ or as a standalone device for your existing system.



Small
27" - 37" Circumference



Large
33" - 51" Circumference



Extra Large
39" - 59" Circumference

COMP-KIT-A	FreedomX™ Compression Belt Kit Includes: Small and Large Belts with Sphygmomanometer	1 ea.
COMP-KIT-B	FreedomX™ Compression Belt Kit Include: Large and Extra-Large Belts with Sphygmomanometer	1 ea.
COMP-KIT-C	FreedomX™ Compression Belt Kit Includes: Small, Large, Extra-Large Belts and Sphygmomanometer	1 ea.
COMP-SPG	FreedomX™ Compression Belt Sphygmomanometer	1 ea.

See cdrsys.ca for more information on the FreedomX Modules and FreedomX Accessories.