




# FREEDOMX™

THE FREEDOMX™ INTUITION CLICK™ MODULE



The Intuition Click™ offers superior comfort and accuracy. Ease and accuracy with you in mind.


## Precision, Non-invasive and Frameless

The Intuition Click™ provides a system to produce a fully customized, rigid, air equivalent, non-invasive positioning solution that provides full posterior support to a patient's head, C-spine, and shoulders. The completed customized support provides rigidity right up to the coronal plane.

Combined with either a closed mask or open-face mask, the Intuition Click™ has been shown to better the performance of the most accurate SRS immobilization solutions with superior patient comfort.



Intuition Click Module



The Intuition Click Module is one of the ten highly-designed FreedomX™ Modules.

The Instaform™ Kit  
Easy to mix, the Instaform™ solutions are quickly combined to create one solution for accurate molding.



See [cdrsys.ca](http://cdrsys.ca) for more information on the FreedomX Modules and FreedomX Accessories.




# FREEDOMX™

THE FREEDOMX™ INTUITION CLICK™ MODULE

## The Advantage of an Intuition Click™ created Mold



The Intuition Mold™ creates a true rigid active hold at the cranial stop.



### 1/ A TRUE MATCH


The Intuition Click™ creates an Intuition Mold™ that is a *true* match to the patient's contour.

### 2/ AN ACTIVE HOLD AT THE CORONAL PLANE

The Intuition Mold™ is rigid up to the coronal plane for customized support that never relaxes. It creates a true active hold with pressure at the axial part to prevent roll and pitch.

### 3/ C-SPINE HOLD

The Intuition Mold™ is created using Instaform™. It forms perfectly around the patients head and shoulders – the result is a C-spine hold that tapers to zero.




## The Ultimate in Accuracy > The Intuition™ Click is a frame used in the patient setup phase to create a custom mold unlike any other!

The Intuition Click Module™ allows you to create a custom Intuition Mold™ that provides full rigid support to a patient's head, neck and shoulders, resulting in superior positioning accuracy.

*The creation of the Intuition Click™ Mold using the Instaform™ can be a transformative experience for the patient. Its warmth calms the patient. It has been described as "feeling like a warm bath."*

The Intuition Mold™ is as light as air and as solid as the ground we walk on.

<b>FDM-S-SRSHN</b>	<b>Intuition Click™ Module for SRS/SRT + Head &amp; Neck Setup &amp; Simulation</b>	1 pkg.
	Includes: MM14, MM14-IR, LB-HP, STA-FILL-X, STA-12-X, THS03-1-MM	
<b>FDM-S-SRS</b>	<b>Intuition Click™ Module for SRS/SRT &amp; Simulation</b>	1 pkg.
	Includes: MM14-IF, MM14-IR, STA-FILL-X, THS03-1-MM	



<b>MM14</b>	Intuition Click™ for Custom Head, Neck & Shoulder includes frame + risers	1 ea.
<b>MM14-FI</b>	Intuition Click™ Frame Replacement Foam Insert Kit	1 ea.
<b>MM14-IR</b>	Intuition Click™ Risers – for use with F size tall head supports	1 ea.
<b>STA-12-X</b>	Shoulder-Tech – compatible with LB-SBRT-X, LB02 & SB02	1 ea.
<b>STA-FILL-X</b>	FreedomX™ Cranial SRS Adapters for use with SB02-S1	1 ea.
<b>LB-HP</b>	FreedomX™ Arms Down Hand Indexer – configured to couch	1 ea.
<b>AAS03</b>	Arm/Elbow Adjustable Support Strap	1 ea.