




# FREEDOMX™

THE FREEDOMX™ MASK STABILIZATION MODULE



Mask Stabilization offers extra stabilization. Designed with you in mind.


## Unparalleled Stabilization

The Mask Stabilization System is an optional add-on that provides additional rigidity and stabilization to meet immobilization protocols. Often used when optimal posterior support is not the focus. This system attaches directly to an immobilization mask and eliminates the need for an invasive bite block system as a means to enhance patient immobilization.

Both the *Wishbone Stabilizer* and *Intuition Stabilizer* index directly to all FreedomX™ Overlays and Couch Extensions using our proprietary LT-Thermoplastics Quick-Attachment process. Both open and closed face masks can be used with FreedomX™ Mask Stabilization System.



### Mask Stabilization Module



The Mask Stabilization Module is one of the ten highly-designed FreedomX™ Modules.

NWB02-MM	Wishbone Carbon Fiber Stabilizer – large – for use with Intuition	1 ea.
NWB02	Wishbone Mask Stabilizer System – for Freedom Overlays	1 ea.
NWB02-QA	Quick Attachments x 12 ea.	1 pkg.

IFF03	Instaform™ Kit – Standard (4 Kits/Box)	4 ea.
IFF03-B	Instaform™ Kit – Large (4 Kits/Box)	4 ea.
IFF03-D	Instaform™ Kit – Extra Large (8 Kits/Box)	8 ea.
IFF03-E	Instaform™ Kit – Treble Size (6 Kits/Box)	6 ea.

See [cdrsys.ca](http://cdrsys.ca) for more information on the FreedomX Modules and FreedomX Accessories.