




# HEXSPAND™

HONEYCOMB LT-THERMOPLASTIC™  
IMMOBILIZATION MASKS + SHEETS

CDR Systems LT-Thermoplastics™ are manufactured using the most current and advanced processes available, ensuring the highest consistency with respect to stretch, working time, and ultimate patient comfort.

With a honeycomb pattern, HEXSPAND™ masks are stronger than before without compromising comfort or functionality.



Hexspand – detail




Shown here:  
< Hexspand Closed Mask  
(blue) – formed example



Hexspand Closed Mask >  
(white) – formed example



Shown here:  
Hexspand Open Mask >  
(blue) – formed example



Shown here:  
< Hexspand Closed Mask  
(white) – flat  
Hexspand Closed Mask >  
(blue) – flat





See [cdrsys.ca](http://cdrsys.ca) for more information on the many HEXSPAND Thermoplastic options available.



**BLUE VERSUS WHITE – THE DIFFERENCE:**

Blue LT-Thermoplastics represent the latest innovation in our immobilization thermoplastic materials. The material blend has been optimized for just the right amount of increased working time and a lower force moulding experience for both user and patient when compared to our White LT-Thermoplastic material.

*The result is increased patient comfort while maintaining the accuracy and rigidity you expect from CDR Systems.*



**LT-THERMOPLASTIC™ MASKS + THE NEW FREEDOM DUO™ SRS IMMOBILIZATION SYSTEM WITH FITADJUST™:**

Combining Thermoplastic masks with a system that provides the ultimate in patient comfort and adjustability that has never been this easy.



LT-Thermoplastic™ Duo Hexspand Mask Set with Shoulder Tech

LT-Thermoplastic™ Duo Hexspand Mask Set



# THERMOPLASTICS

## LT-THERMOPLASTIC™ DRAPE + MOULD MASKS IMMOBILIZATION MASKS

CDR Systems LT-Thermoplastics™ are manufactured using the most current and advanced processes available, ensuring the highest consistency with respect to stretch, working time, and ultimate patient comfort.


With a honeycomb pattern, HEXSPAND™ masks are stronger than before without compromising comfort or functionality.  
**Designed with you in mind!**




LT-Thermoplastic™  
Blue Hexspand Mask




LT-Thermoplastic™  
White Hexspand Open Mask



Drape + Mould LT-Thermoplastic™  
(Shown here: Hexspand Extra Long  
H24XLX/H32XLX)




Open Drape + Mould LT-Thermoplastic™  
(Shown here: Standard  
H240/H320)



Open Drape + Mould LT-Thermoplastic™  
(Shown here: Hexspand Hybrid  
H24Y0/H32Y0)

\*Hybrid Open Face Mask: Engineered for ultimate patient comfort. This mask never covers the patient's face, making it especially ideal for claustrophobic patients.



Drape + Mould LT-Thermoplastic™  
(Shown here: Reinforced  
H24RO/H32RO)

2.4 MM

<b>H24</b>	Drape and Mould LT-Thermoplastic™ – White Standard	5 pk.
<b>H24X</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>H24BL</b>	Drape and Mould LT-Thermoplastic™ – Blue Standard	5 pk.
<b>H24BLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>H24L</b>	Drape and Mould LT-Thermoplastic™ – White Long	5 pk.
<b>H24LX</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>H24LBL</b>	Drape and Mould LT-Thermoplastic™ – Blue Long	5 pk.
<b>H24LBLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>H24XL</b>	Drape and Mould LT-Thermoplastic™ – White Extra Long	5 pk.
<b>H24XLX</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Extra Long	5 pk.
<b>H24XLBL</b>	Drape and Mould LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>H24XLBLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.
<b>H240</b>		
<b>H240X</b>	Open Drape and Mould LT-Thermoplastic™ – White Standard	5 pk.
<b>H240BL</b>	Open Drape and Mould LT-Thermoplastic™ – White Hexspand	5 pk.
<b>H240BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Standard	5 pk.
<b>H240BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>H24Y0</b>	Open Drape and Mould LT-Thermoplastic™ – White Hybrid	5 pk.
<b>H24Y0X</b>	Open Drape and Mould LT-Thermoplastic™ – White Hexspand Hybrid	5 pk.
<b>H24Y0BL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hybrid	5 pk.
<b>H24Y0BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hexspand Hybrid	5 pk.
<b>H24Y0BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Black Hexspand Hybrid	5 pk.
<b>H24RO</b>	Open Drape and Mould LT-Thermoplastic™ – White Reinforced	5 pk.
<b>H24ROBL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Reinforced	5 pk.
<b>H24ROBK</b>	Open Drape and Mould LT-Thermoplastic™ – Black Reinforced	5 pk.

3.2 MM

<b>H32</b>	Drape and Mould LT-Thermoplastic™ – White Standard	5 pk.
<b>H32X</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>H32BL</b>	Drape and Mould LT-Thermoplastic™ – Blue Standard	5 pk.
<b>H32BLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>H32L</b>	Drape and Mould LT-Thermoplastic™ – White Long	5 pk.
<b>H32LX</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>H32LBL</b>	Drape and Mould LT-Thermoplastic™ – Blue Long	5 pk.
<b>H32LBLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>H32XL</b>	Drape and Mould LT-Thermoplastic™ – White Extra Long	5 pk.
<b>H32XLX</b>	Drape and Mould LT-Thermoplastic™ – White Hexspand Extra Long	5 pk.
<b>H32XLBL</b>	Drape and Mould LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>H32XLBLX</b>	Drape and Mould LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.
<b>H320</b>		
<b>H320X</b>	Open Drape and Mould LT-Thermoplastic™ – White Standard	5 pk.
<b>H320X</b>	Open Drape and Mould LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>H320BL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Standard	5 pk.
<b>H320BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>H32Y</b>	Open Drape and Mould LT-Thermoplastic™ – White Hybrid	5 pk.
<b>H32YX</b>	Open Drape and Mould LT-Thermoplastic™ – White Hexspand Hybrid	5 pk.
<b>H32YBL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hybrid	5 pk.
<b>H32YBLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hexspand Hybrid	5 pk.
<b>H32Y0</b>	Open Drape and Mould LT-Thermoplastic™ – White Hybrid	5 pk.
<b>H32Y0X</b>	Open Drape and Mould LT-Thermoplastic™ – White Hexspand Hybrid	5 pk.
<b>H32Y0BL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hybrid	5 pk.
<b>H32Y0BLX</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Hexspand Hybrid	5 pk.
<b>H32RO</b>	Open Drape and Mould LT-Thermoplastic™ – White Reinforced	5 pk.
<b>H32ROBL</b>	Open Drape and Mould LT-Thermoplastic™ – Blue Reinforced	5 pk.



# THERMOPLASTICS

## LT-THERMOPLASTIC™ HEAD, NECK & SHOULDER MASKS IMMOBILIZATION MASKS




Head, Neck & Shoulder LT-Thermoplastic™  
Blue Hexspand Open

### HEAD, NECK & SHOULDER LT-THERMOPLASTIC™ > IMMOBILIZATION MASKS

#### 2.4 MM

<b>N24</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Standard	5 pk.
<b>N24X</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>N24BL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Standard	5 pk.
<b>N24BLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>N240</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Open	5 pk.
<b>N240X</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Open	5 pk.
<b>N240BL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Open	5 pk.
<b>N240BLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Open	5 pk.
<b>N24YO</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hybrid Open	5 pk.
<b>N24YOBL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hybrid Open	5 pk.
<b>N24PEX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Pediatric	5 pk.
<b>N24PEBLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Pediatric	5 pk.




Head, Neck & Shoulder LT-Thermoplastic™  
(Shown here: Open N240/N320)


#### 3.2 MM

<b>N32</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Standard	5 pk.
<b>N32X</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>N32BL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Standard	5 pk.
<b>N32BLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>N320</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Open	5 pk.
<b>N320X</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Open	5 pk.
<b>N320BL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Open	5 pk.
<b>N320BLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Open	5 pk.
<b>N32YO</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hybrid Open	5 pk.
<b>N32YOBL</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hybrid Open	5 pk.
<b>N32PEX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – White Hexspand Pediatric	5 pk.
<b>N32PEBLX</b>	Head, Neck & Shoulder LT-Thermoplastic™ – Blue Hexspand Pediatric	5 pk.


**SLIM U-FRAME LT-THERMOPLASTIC™ > IMMOBILIZATION MASKS**



Slim U-Frame LT-Thermoplastic™  
(Shown here: Hexspand Long  
S24LX/S32LX)



Slim U-Frame LT-Thermoplastic™  
(Shown here: Extra Long  
S24/S32)



Split U-Frame LT-Thermoplastic™  
(Shown here: Long  
SP24X/SP32X)

**2.4 MM**

<b>S24</b>	Slim U-Frame LT-Thermoplastic™ – White Standard	5 pk.
<b>S24X</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>S24BL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Standard	5 pk.
<b>S24BLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>S24L</b>	Slim U-Frame LT-Thermoplastic™ – White Long	5 pk.
<b>S24LX</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>S24LBL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Long	5 pk.
<b>S24LBLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>S24XL</b>	Slim U-Frame LT-Thermoplastic™ – White Extra Long	5 pk.
<b>S24XLX</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Extra Long	5 pk.
<b>S24XLBL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>S24XLBLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.

**3.2 MM**

<b>S32</b>	Slim U-Frame LT-Thermoplastic™ – White Standard	5 pk.
<b>S32X</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>S32BL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Standard	5 pk.
<b>S32BLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>S32L</b>	Slim U-Frame LT-Thermoplastic™ – White Long	5 pk.
<b>S32LX</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>S32LBL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Long	5 pk.
<b>S32LBLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>S32XL</b>	Slim U-Frame LT-Thermoplastic™ – White Extra Long	5 pk.
<b>S32XLX</b>	Slim U-Frame LT-Thermoplastic™ – White Hexspand Extra Long	5 pk.
<b>S32XLBL</b>	Slim U-Frame LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>S32XLBLX</b>	Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.

**SPLIT SLIM U-FRAME LT-THERMOPLASTIC™ > IMMOBILIZATION MASKS**

**2.4 MM**

<b>SP24</b>	Split Slim U-Frame LT-Thermoplastic™ – Standard	5 pk.
<b>SP24X</b>	Split Slim U-Frame LT-Thermoplastic™ – Hexspand Standard	5 pk.
<b>SP24BL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Standard	5 pk.
<b>SP24BLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>SP24L</b>	Split Slim U-Frame LT-Thermoplastic™ – White Long	5 pk.
<b>SP24LX</b>	Split Slim U-Frame LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>SP24LBL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Long	5 pk.
<b>SP24LBLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>SP24XL</b>	Split Slim U-Frame LT-Thermoplastic™ – Extra Long	5 pk.
<b>SP24XLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Hexspand Extra Long	5 pk.
<b>SP24XLBL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>SP24XLBLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.

**3.2 MM**

<b>SP32</b>	Split Slim U-Frame LT-Thermoplastic™ – White Standard	5 pk.
<b>SP32X</b>	Split Slim U-Frame LT-Thermoplastic™ – White Hexspand Standard	5 pk.
<b>SP32BL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Standard	5 pk.
<b>SP32BLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Standard	5 pk.
<b>PS32L</b>	Split Slim U-Frame LT-Thermoplastic™ – White Long	5 pk.
<b>SP32LX</b>	Split Slim U-Frame LT-Thermoplastic™ – White Hexspand Long	5 pk.
<b>SP32LBL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Long	5 pk.
<b>SP32LBLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Long	5 pk.
<b>SP32XL</b>	Split Slim U-Frame LT-Thermoplastic™ – White Extra Long	5 pk.
<b>SP32XLX</b>	Split Slim U-Frame LT-Thermoplastic™ – White Hexspand Extra Long	5 pk.
<b>SP32XLBL</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Extra Long	5 pk.
<b>SP32XLBLX</b>	Split Slim U-Frame LT-Thermoplastic™ – Blue Hexspand Extra Long	5 pk.

LT-Thermoplastic™ perforated  
12" x 24" chest mold >





## SABELLA FLEX BREAST LT-THERMOPLASTIC™

### 2.4 MM


<b>BS24P</b>	Sabella Flex Breast LT-Thermoplastic™ – White Perforated	1 ea.
<b>BS24-01</b>	Sabella Flex Breast LT-Thermoplastic™ – White Solid	1 ea.
<b>BS24X</b>	Sabella Flex Breast LT-Thermoplastic™ – White Hexspand	1 ea.
<b>BS24BLP</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Perforated	1 ea.
<b>BS24BL-01</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Solid	1 ea.
<b>BS24BLX</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Hexspand	1 ea.

### 3.2 MM

<b>BS32P</b>	Sabella Flex Breast LT-Thermoplastic™ – White Perforated	1 ea.
<b>BS32-01</b>	Sabella Flex Breast LT-Thermoplastic™ – White Solid	1 ea.
<b>BS32X</b>	Sabella Flex Breast LT-Thermoplastic™ – White Hexspand	1 ea.
<b>BS32BLP</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Perforated	1 ea.
<b>BS32BL-01</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Solid	1 ea.
<b>BS32BLX</b>	Sabella Flex Breast LT-Thermoplastic™ – Blue Hexspand	1 ea.



Hexspand Closed Head Mask  
with White Hexspand Shoulder Tech >




## SHOULDER TECH LT-THERMOPLASTIC™

### 2.4 MM

<b>ST24-01</b>	Shoulder Tech LT-Thermoplastic™ – White	5 pk.
<b>ST24X</b>	Shoulder Tech LT-Thermoplastic™ – White Hexspand	5 pk.
<b>ST24BL-01</b>	Shoulder Tech LT-Thermoplastic™ – Blue	5 pk.
<b>ST24BLX</b>	Shoulder Tech LT-Thermoplastic™ – Blue Hexspand	5 pk.

### 3.2 MM

<b>ST32-01</b>	Shoulder Tech LT-Thermoplastic™ – White	5 pk.
<b>ST32X</b>	Shoulder Tech LT-Thermoplastic™ – White Hexspand	5 pk.
<b>ST32BL-01</b>	Shoulder Tech LT-Thermoplastic™ – Blue	5 pk.
<b>ST32BLX</b>	Shoulder Tech LT-Thermoplastic™ – Blue Hexspand	5 pk.



Shoulder Tech LT-Thermoplastic™  
(Shown here: ST32-01)



# THERMOPLASTICS

## LT-THERMOPLASTIC™ PERFORATED AND SOLID SHEETS

All LT-Thermoplastic sheets are manufactured using the most advanced processes ensuring the highest consistency with respect to stretch, working time and ultimate patient comfort. The sheets come in various sizes and thicknesses and are available in both Solid and Perforated types to ensure maximum accuracy and comfort.

Simply attach a set of reusable handles to the LT-Thermoplastic sheet and attach to the frame (shown here on the SaBella Flex™). Once treatment is complete, retain the handles for reuse and recycle the used LT-Thermoplastic™.



### TECHNICAL SPECIFICATIONS

#### ACTIVATION TEMPERATURE

Between 70°C and 75°C  
(160°F to 165°F)



#### HEATING TIME

Water Bath: 1-2 minutes  
(max 5 min)

Dry Heat Oven: 3-15 minutes  
(max 30 min)

MOLDING TIME: 1-2 minutes

HARDENING TIME: min 10 minutes



All LT-Thermoplastic products – including Perforated and Solid Sheets – are dry heat compatible.

## FREEDOM PELVIS & THORAX LT-THERMOPLASTIC™

### 3.2 MM

<b>PLV-1232P</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Perforated, 12" x 18"	1 ea.
<b>PLV-1232-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Solid, 12" x 18"	1 ea.
<b>PLV-1232BLP</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Perforated, 12" x 18"	1 ea.
<b>PLV-1232BL-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Solid, 12" x 18"	1 ea.
<b>PLV-1832P</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Perforated, 18" x 18"	1 ea.
<b>PLV-1832-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Solid, 18" x 18"	1 ea.
<b>PLV-1832BLP</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Perforated, 18" x 18"	1 ea.
<b>PLV-1832BL-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Solid, 18" x 18"	1 ea.
<b>PLV-2432P</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Perforated, 24" x 18"	1 ea.
<b>PLV-2432-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – White Solid, 24" x 18"	1 ea.
<b>PLV-2432BLP</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Perforated, 24" x 18"	1 ea.
<b>PLV-2432BL-01</b>	Freedom Pelvis & Thorax LT-Thermoplastic™ – Blue Solid, 24" x 18"	1 ea.

## SHEET LT-THERMOPLASTIC™

### 1.6 MM

<b>RS16-12x18</b>	Sheet LT-Thermoplastic™ – White Perforated, 12" x 18"	1 ea.
<b>RS16-12x18-01</b>	Sheet LT-Thermoplastic™ – White Solid, 12" x 18"	1 ea.
<b>RS16-12x18BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 12" x 18"	1 ea.
<b>RS16-12x18BL-01</b>	Sheet LT-Thermoplastic™ – Blue Solid, 12" x 18"	1 ea.
<b>RS16-18x24</b>	Sheet LT-Thermoplastic™ – White Perforated, 18" x 24"	1 ea.
<b>RS16-18x24-01</b>	Sheet LT-Thermoplastic™ – White Solid, 18" x 24"	1 ea.
<b>RS16-18x24BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 18" x 24"	1 ea.
<b>RS16-18x24BL-01</b>	Sheet LT-Thermoplastic™ – Blue Solid, 18" x 24"	1 ea.

### 2.0 MM


<b>RS20-12x18</b>	Sheet LT-Thermoplastic™ – White Perforated, 12" x 18"	1 ea.
<b>RS20-12x18BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 12" x 18"	1 ea.
<b>RS20-18x24</b>	Sheet LT-Thermoplastic™ – White Perforated, 18" x 24"	1 ea.
<b>RS20-18x24-01</b>	Sheet LT-Thermoplastic™ – White Solid, 18" x 24"	1 ea.
<b>RS20-18x24BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 18" x 24"	1 ea.
<b>RS20-18x24BL-01</b>	Sheet LT-Thermoplastic™ – Blue Solid, 18" x 24"	1 ea.

### 3.2 MM

<b>RS32-12x18</b>	Sheet LT-Thermoplastic™ – White Perforated, 12" x 18"	1 ea.
<b>RS32-12x18BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 12" x 18"	1 ea.
<b>RS32-18x24</b>	Sheet LT-Thermoplastic™ – White Perforated, 18" x 24"	1 ea.
<b>RS32-18x24-01</b>	Sheet LT-Thermoplastic™ – White Solid, 18" x 24"	1 ea.
<b>RS20-18x24BL</b>	Sheet LT-Thermoplastic™ – Blue Perforated, 18" x 24"	1 ea.
<b>RS20-18x24BL-01</b>	Sheet LT-Thermoplastic™ – Blue Solid, 18" x 24"	1 ea.




## S-TYPE MASKS



ST-STYP-24 + 32

ST-STYP-24R + 32R

<b>ST-STYP-24</b>	LT-Thermoplastic S Type Immobilization Mask Head & Shoulder – 2.4mm Standard	5 ea.
<b>ST-STYP-32</b>	LT-Thermoplastic S Type Immobilization Mask Head & Shoulder – 3.2mm Standard	5 ea.
<b>ST-STYP-24R</b>	LT-Thermoplastic S Type Immobilization Mask Head & Shoulder – 2.4mm Reinforced	5 ea.
<b>ST-STYP-32R</b>	LT-Thermoplastic S Type Immobilization Mask Head & Shoulder – 3.2mm Reinforced	5 ea.




SH-STYP-24 + 32

SH-STYP-24R + 32R

<b>SH-STYP-24</b>	LT-Thermoplastic S Type Immobilization Mask Head Only – 2.4mm Standard	5 ea.
<b>SH-STYP-32</b>	LT-Thermoplastic S Type Immobilization Mask Head Only – 3.2mm Standard	5 ea.
<b>SH-STYP-24R</b>	LT-Thermoplastic S Type Immobilization Mask Head Only – 2.4mm Reinforced	5 ea.
<b>SH-STYP-32R</b>	LT-Thermoplastic S Type Immobilization Mask Head Only – 3.2mm Reinforced	5 ea.

## BL TYPE IMMOBILIZATION MASKS NEW IMAGES



BL410


BL408

BL406

<b>BL410</b>	LT-Thermoplastic BL-Type Immobilization Mask Mask Head, Neck & Shoulder Type	5 ea.
<b>BL408</b>	LT-Thermoplastic BL-Type Immobilization Mask Mask B-Type	5 ea.
<b>BL406</b>	LT-Thermoplastic BL-Type Immobilization Mask Mask A-Type	5 ea.

## BATH LINERS + TRANSFER SHEETS + SHIMS

### Water Bath Liners & Transfer Sheets



<b>99-01-S</b>	Small liners for water bath and dry heat applications (non-stick) 10" x 15" / 254mm x 381mm	10 ea.
<b>99-01-M</b>	Medium liners for water and dry heat application (non-stick) 18" x 24" / 457mm x 609mm	10 ea.
<b>99-01-L</b>	Larger liners for water bath and dry heat applications (non-stick) 24" x 36" / 609mm x 914mm	10 ea.

### Shims



<b>MSHIM</b>	Face Mask Shims Thickness: 0.5mm	6 ea.
<b>MSHIM-I</b>	Intuition Click™ Shims Thickness: 0.5mm	3 ea.
<b>HSIM</b>	Silverman Head Support Shims Thickness: 0.5mm	3 ea.