



# VACUUM LOCK™

## CUSHIONS

The Vacuum Lock™ Cushions are structurally reinforced nylon with double welded seam giving ultimate strength and durability.

**ADDITIONAL FEATURES:**

- > Stepped valve for secure vacuum hold and ease of use
- > Ultra low density bead fill provides optimal patient conformity
- > Proven and tested to hold shape through course of treatment.
- > Our innovative High Friction Hybrid option provides amplified indexing and locking on any surface.



STANDARD

HIGH FRICTION HYBRID



Vacuum Lock™ Cushion shown on FreedomX™ Total Body Precision Patient Immobilization System, Arms Up Configuration (LB-SBRT).

See [cdrsys.ca](http://cdrsys.ca) for more information on the Vacuum Lock.

VB-4080F15   VBF-4080F15	Vacuum Lock™ Cushion – 40 x 80cm 15L fill for extremities	1 ea.
VB-5070F15   VBF-5070F15	Vacuum Lock™ Cushion – 50 x 70cm 15L fill for extremities	1 ea.
VB-7070F20   VBF-7070F20	Vacuum Lock™ Cushion – 70 x 70cm 20L fill for upper body	1 ea.
VB-70150F60   VBF-70150F60	Vacuum Lock™ Cushion – 70 x 150cm 60L fill for hip or pelvis	1 ea.
VB-80150F65   VBF-80150F65	Vacuum Lock™ Cushion – 80 x 150cm 65L fill for hip or pelvis	1 ea.
VB-10070F30   VBF-10070F30	Vacuum Lock™ Cushion – 100 x 70cm 30L fill for hip, pelvis & upper body	1 ea.
VB-10070F35   VBF-10070F35	Vacuum Lock™ Cushion – 100 x 70cm 35L fill for hip, pelvis & upper body	1 ea.
VB-100100F30   VBF-100100F30	Vacuum Lock™ Cushion – 100 x 100cm 30L fill for upper body	1 ea.
VB-100150F70   VBF-100150F70	Vacuum Lock™ Cushion – 100 x 150cm 70L fill for pelvis or full body	1 ea.
VB-100220F120   VBF-100220F120	Vacuum Lock™ Large – 100 x 220cm 120L fill for whole body	1 ea.
VB-T10070X2F30   VBF-T10070X2F30	Vacuum Lock™ T Shape – Dual Chamber T configuration (2 chambers 70 x 100cm 30L fill)	1 ea.
VB-T6565F7   VBF-T6565F7	Vacuum Lock™ T Shape – 65 x 65cm 7L fill T shaped for head and shoulders	1 ea.
VB-100130PENF45   VBF-100130PENF45	Vacuum Lock™ Pentagon – 100 x 130cm 45L fill pentagon shape	1 ea.
VB-T7993F20   VBF-T7993F20	Vacuum Lock™ T Shape – 79 x 93cm 20L fill for SaBella Flex™	1 ea.
VB-20-014	Vacuum Lock™ Cover – Valve Protector for Vacuum Lock™ Cushions	1 ea.


VB = Standard    VBF = High Friction Hybrid



# VACUUM LOCK™

## INDEX BARS + ACCESSORIES

A single Indexer may be used on any one cushion. An additional indexer may be used as required for increased reposition accuracy on longer cushion.



North America Only



<b>IB06-BA</b>	<b>Twin Index System</b> for vacuum bags. Attaches to a standard indexing bar	1 ea.
<b>IB06-BA-BX</b>	<b>Quad Index System</b> for vacuum bags. Attaches to a standard indexing bar	1 ea.

<b>VB-HOSE</b>	Vacuum Lock™ cushion adapter hose	1 ea.
<b>VB-HOSE-G</b>	Vacuum Lock™ cushion hose with pressure gauge	1 ea.

<b>VBPUMP3</b>	<b>Vacuum Lock™ 100/120V</b> – heavy duty vacuum pump with manual inflate/deflate connections and built in vacuum gauge. Includes hose attachment ( <b>VB-HOSE</b> ).	1 ea.
<b>VBPUMP3-EU</b>	<b>Vacuum Lock™ 220/240V</b> – heavy duty vacuum pump with manual inflate/deflate connections and built in vacuum gauge. Includes hose attachment ( <b>VB-HOSE</b> ).	1 ea.

<b>VBPUMP2</b>	<b>Bench or cart-based Vacuum Lock™ pump – 100/120V</b> , foot pedal operated inflate/deflate with audible pressure alarm. Includes hose attachment ( <b>VB-HOSE</b> ).	1 ea.
----------------	---	-------

<b>VBPUMP1</b>	<b>Portable vacuum pump</b> – manually adjusted inflate/deflate connections. Includes hose attachment ( <b>VB-HOSE</b> )	1 ea.
<b>VBPUMP1-EU</b>	<b>Portable vacuum pump 220/240 Volt compatible</b> . Includes hose attachment ( <b>VB-HOSE</b> ).	1 ea.

<b>VBMSC</b>	<b>Vacuum Lock™ Mobile Hanging Rail Cart</b> – Stores multiple vacuum cushions safely and rolls on wheels, adjustable height in three positions, 140 cm, 152 cm and 165 cm.	1 ea.
<b>VBHOOKS</b>	<b>Vacuum Lock™ Hooks</b> – 10 hooks with rubber protected ends	1 pkg.