



《EXTREMITY + TOTIM》

Dedicated Extremity Immobilization and Posterior Support System
 Your New Extremity Immobilization Solution



Low Shrinkage




SGRT Compatible



Non-Stick



Latex Free



Low activation temperature



X-Ray translucent



Phthalates Free




MR Compatible available

CDR's new solution for Patient Immobilization, designed specifically for Extremity Indexing which ensures secure, precise positioning for challenging anatomical areas.

Engineered for superior rigidity and reproducibility, the new TOTIM™ Customized Patient Posterior Cushions support and contour to the patients unique shape and is then indexed using a LT-Thermoplastic Hexspand™ mask. This directly addresses the common difficulties of efficient reproducibility and achieving a comfortable patient experience.

Extremity Immobilization w. Customized Posterior Positioning



Shown here:
 ▲ **EB24-STD**
 Extremity Board for E-Type Thermoplastics



Shown here:
 ▲ **TTM3010XS** (20x25cm)
 TOTIM Custom Posterior Cushion - XS




Shown here:
 ▲ **TTM3010S** (25x35cm)
 TOTIM Custom Posterior Cushion - S




Shown here:
 ▲ **TTM3006** (70x30cm)
 TOTIM Custom Posterior Cushion - Limb




Shown here:
 ▲ **E-P32X3P**
 E-Type Hexspand 3 Point Mask
 (Best for Foot and Ankle Immobilization)



Shown here:
 ▲ **E-P32X4P**
 E-Type Hexspand 4 Point Mask
 (Best for Hand, Knee, and Leg)



Shown here:
 ▲ **BS32P**
 E-Type Hexspand 4-point Mask DH 3.2mm
 (Best for Double Hand)



Shown here:
 ▲ **E-P32X4PE**
 E-Type Hexspand 4 Point Elbow Mask
 (Best for Elbow or Arm Immobilization)



Extremity Immobilization w. Customized Posterior Positioning



Shown here ^ E-P32X3P & TTM3006




Shown here ^ E-P32X4P & TTM3010XS



Shown here ^ E-P32X4PE & TTM3006



Shown here ^ E-P32X3P, TTM3006, EB-24STD



Shown here ^ BS32PX, EB-24STD



Shown here ^ E-P32X4PE, TTM3006, EB-24STD

Packages
EXT-PKG

Comprehensive Package for all Extremity Immobilization
Includes: 1x EB-24STD, 1x E-P32X4P, 1x E-P32X4PE, 1x E-P32X3P, 1x BS32P, 1x TTM3010XS, 1x TTM3006, 1x E-PFCC, 1x E-PIK

PKG

Extremity

EB-24STD
E-P32X4P
E-P32X4PE
E-P32X3P
BS32PX

Extremity Overlay for E-Type Thermoplastics
E-Type Extremity Hexspand 4 Point Mask 3.2mm
E-Type Extremity Hexspand 4 Point Elbow Mask 3.2mm
E-Type Extremity Hexspand 3 Point Mask 3.2mm
LT-Thermoplastic Hexspand 4-point Gen II Mask - Dual Hand 3.2mm

EA
5 PK
5 PK
5 PK
5 PK

Positioning Aids

TTM3010XS
TTM3010S
TTM3006
E-PFCC
E-PIK

TOTIM Customized Patient Posterior Cushion - XS
TOTIM Customized Patient Posterior Cushion - S
TOTIM Customized Patient Posterior Cushion - LIMB
E-Type Extremity Foam Comfort Cushions Kit (Contains 5x Small, 5x Large EVA Pads)
E-Type Extremity Indexing Kit for TOTIM Cushions (Contains 3x Indexing Plate)

10 PK
10 PK
2 PK
10 PK
3 PK

Contact your sales rep more information or to schedule a demo. You can also direct inquiries to info@cdrsys.ca